




# *VENA CRYL 910*



**vena**

Vena Medical Suppliers  
[www.sutures.com.eg](http://www.sutures.com.eg)

# VENA CRYL 910

ABSORBABLE SURGICAL SUTURE  
POLYGLACTINE 910 BRAIDED



## Type

Braided and coated multifilament

## Composition

Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide

## Coating

Poly-glycolide-co-L-lactide and calcium stearate (<1%)

## Colour

Violet

## Tissue reaction

Minimal

## Absorption

The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days  
Between 40% and 50% of tensile strength remain after 21 days.

## BENEFITS

*Excellent tissue support*

*Minimal tissue reaction*

*Predictable absorption rate*

*Excellent handling properties*

*High flexibility*

*High tensile strength*

*Excellent knot security*

*Soft passage through the tissues  
with minimal sawing,  
tissue drag and trauma*

*Good knot tie down / adjustability*

*Reduced capillary action*

*Guaranteed seal and product sterility*

## CHARACTERISTICS

*High tensile strength*

*Synthetic*

*Breaks down by hydrolysis*

*Braided structure*

*Coated multifilament*

*Hermetically sealed packing*

## INDICATIONS

*General soft tissue closing  
and/or ligation*

*General surgery*

*Skin closure / Gastrointestinal surgery*


*Plastic surgery / Ophthalmic surgery*

*Gynaecology / Obstetrics*

*Orthopaedics / Urology*

 **VENA CRYL 910**

ABSORBABLE SURGICAL SUTURE  
POLYGLACTINE 910 BRAIDED




ROUND BODIED TAPER POINT - 1/2 CIRCLE

| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>OHR 17 mm</b><br>   | 1   | 5/0 | 75 cm  | violet | VC50OHR17    |
|                                                                                                         | 1.5 | 4/0 | 75 cm  | violet | VC40OHR17    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30OHR17    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHR17    |
| <b>OHR 22 mm</b><br>  | 1.5 | 4/0 | 75 cm  | violet | VC40OHR22    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30OHR22    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHR22    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01OHR22    |
| <b>OHR 26 mm</b><br> | 1.5 | 4/0 | 75 cm  | violet | VC40OHR26    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30OHR26    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHR26    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00OHR26    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01OHR26    |
|                                                                                                         | 5   | 2   | 75 cm  | violet | VC02OHR26    |
| <b>OHR 30 mm</b><br> | 2   | 3/0 | 75 cm  | violet | VC30OHR30    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHR30    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00OHR30    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01OHR30    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHR30    |
| <b>OHR 36 mm</b><br> | 2   | 3/0 | 75 cm  | violet | VC30OHR36    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHR36    |
|                                                                                                         | 3.5 | 0   | 90 cm  | violet | VC00OHR36    |
|                                                                                                         | 4   | 1   | 90 cm  | violet | VC01OHR36    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHR36    |
| <b>OHR 40 mm</b><br> | 3   | 2/0 | 75 cm  | violet | VC20OHR40    |
|                                                                                                         | 3.5 | 0   | 90 cm  | violet | VC00OHR40    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01OHR40    |
|                                                                                                         | 4   | 1   | 90 cm  | violet | VC01OHR40-90 |
|                                                                                                         | 5   | 2   | 75 cm  | violet | VC02OHR40    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHR40-90 |
| <b>OHR 45 mm</b><br> | 5   | 2   | 90 cm  | violet | VC02OHR45    |

ABSORBABLE SUTURES

 **VENA CRYL 910**

ABSORBABLE SURGICAL SUTURE  
POLYGLACTINE 910 BRAIDED



ROUND BODIED TAPER POINT - HEAVY NEEDLE - 1/2 CIRCLE

| NEEDLE                                                                                                | EP | USP | Length | Colour | Order Number |
|-------------------------------------------------------------------------------------------------------|----|-----|--------|--------|--------------|
| <b>OHR 36 mm</b><br> | 4  | 1   | 90 cm  | violet | VC01OHR36-H  |

ROUND BODIED TAPER POINT - 5/8 CIRCLE

| NEEDLE                                                                                                  | EP | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|----|-----|--------|--------|--------------|
| <b>FER 27 mm</b><br> | 2  | 3/0 | 75 cm  | violet | VC30TER22    |
|                                                                                                         | 3  | 2/0 | 75 cm  | violet | VC20FER27    |


ROUND BODIED TAPER POINT - 3/8 CIRCLE

| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>TER 17 mm</b><br> | 1.5 | 4/0 | 45 cm  | violet | VC40TER17    |
| <b>TER 22 mm</b><br> | 1.5 | 4/0 | 75 cm  | violet | VC40TER22    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30TER22    |
| <b>TER 30 mm</b><br> | 4   | 1   | 75 cm  | violet | VC01TER30    |

REVERSE CUTTING - 1/2 CIRCLE

| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>OHX 17 mm</b><br> | 1   | 5/0 | 75 cm  | violet | VC50OHX17    |
|                                                                                                         | 1.5 | 4/0 | 75 cm  | violet | VC40OHX17    |
| <b>OHX 22 mm</b><br> | 2   | 3/0 | 75 cm  | violet | VC30OHX22    |
| <b>OHX 26 mm</b><br> | 2   | 3/0 | 75 cm  | violet | VC30OHX26    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20OHX26    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00OHX26    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01OHX26    |

ABSORBABLE SUTURES



REVERSE CUTTING - 1/2 CIRCLE


| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>OHX 36 mm</b><br>   | 3   | 2/0 | 75 cm  | violet | VC20OHX36    |
|                                                                                                         | 3.5 | 0   | 90 cm  | violet | VC00OHX36    |
|                                                                                                         | 4   | 1   | 90 cm  | violet | VC01OHX36    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHX36    |
| <b>OHX 40 mm</b><br>  | 2   | 3/0 | 90 cm  | violet | VC30OHX40    |
|                                                                                                         | 3   | 2/0 | 90 cm  | violet | VC20OHX40    |
|                                                                                                         | 3.5 | 0   | 90 cm  | violet | VC00OHX40    |
|                                                                                                         | 4   | 1   | 90 cm  | violet | VC01OHX40    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHX40    |
| <b>OHX 45 mm</b><br> | 4   | 1   | 90 cm  | violet | VC01OHX45    |
|                                                                                                         | 5   | 2   | 90 cm  | violet | VC02OHX45    |

REVERSE CUTTING - 3/8 CIRCLE

| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>TEX 16 mm</b><br> | 1   | 5/0 | 75 cm  | violet | VC50TEX16    |
|                                                                                                         | 1   | 5/0 | 75 cm  | violet | VC50TEX16    |
|                                                                                                         | 1.5 | 4/0 | 75 cm  | violet | VC40TEX16    |
| <b>TEX 19 mm</b><br> | 1   | 5/0 | 75 cm  | violet | VC50TEX19    |
|                                                                                                         | 1.5 | 4/0 | 75 cm  | violet | VC40TEX19    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30TEX19    |
| <b>TEX 24 mm</b><br> | 1.5 | 4/0 | 75 cm  | violet | VC40TEX24    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30TEX24    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20TEX24    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00TEX24    |
| <b>TEX 26 mm</b><br> | 1.5 | 4/0 | 75 cm  | violet | VC40TEX26    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30TEX26    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20TEX26    |
|                                                                                                         | 3   | 2/0 | 90 cm  | violet | VC20TEX26-90 |

 **VENA CRYL 910**

ABSORBABLE SURGICAL SUTURE  
POLYGLACTINE 910 BRAIDED



REVERSE CUTTING - 3/8 CIRCLE

| NEEDLE                                                                                                  | EP  | USP | Length | Colour | Order Number |
|---------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>TEX 26 mm</b><br>   | 1.5 | 4/0 | 75 cm  | violet | VC40TEX26    |
|                                                                                                         | 2   | 3/0 | 75 cm  | violet | VC30TEX26    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20TEX26    |
|                                                                                                         | 3   | 2/0 | 90 cm  | violet | VC20TEX26-90 |
| <b>TEX 30 mm</b><br>   | 2   | 3/0 | 75 cm  | violet | VC30TEX30    |
|                                                                                                         | 3   | 2/0 | 75 cm  | violet | VC20TEX30    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00TEX30    |
| <b>TEX 36 mm</b><br> | 3   | 2/0 | 75 cm  | violet | VC20TEX36    |
|                                                                                                         | 3.5 | 0   | 75 cm  | violet | VC00TEX36    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01TEX36    |
| <b>TEX 40 mm</b><br> | 3.5 | 0   | 75 cm  | violet | VC00TEX40    |
|                                                                                                         | 4   | 1   | 75 cm  | violet | VC01TEX40    |

ROUND BODIED TAPERCUT POINT - 1/2 CIRCLE

| NEEDLE                                                                                                   | EP  | USP | Length | Colour | Order Number |
|----------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>OHRT 26 mm</b><br> | 1.5 | 4/0 | 75 cm  | violet | VC40OHRT26   |
|                                                                                                          | 2   | 3/0 | 75 cm  | violet | VC30OHRT26   |
|                                                                                                          | 3   | 2/0 | 75 cm  | violet | VC20OHRT26   |
| <b>OHRT 30 mm</b><br> | 3   | 2/0 | 75 cm  | violet | VC20OHRT30   |
|                                                                                                          | 3.5 | 0   | 75 cm  | violet | VC00OHRT30   |
|                                                                                                          | 4   | 1   | 75 cm  | violet | VC01OHRT30   |
| <b>OHRT 36 mm</b><br> | 3   | 2/0 | 90 cm  | violet | VC20OHRT36   |
|                                                                                                          | 3.5 | 0   | 90 cm  | violet | VC00OHRT36   |
|                                                                                                          | 4   | 1   | 90 cm  | violet | VC01OHRT36   |
| <b>OHRT 40 mm</b><br> | 3.5 | 0   | 90 cm  | violet | VC00OHRT40   |
|                                                                                                          | 4   | 1   | 90 cm  | violet | VC01OHRT40   |
|                                                                                                          | 5   | 2   | 90 cm  | violet | VC02OHRT40   |


ABSORBABLE SUTURES

 **VENA CRYL 910**

ABSORBABLE SURGICAL SUTURE  
POLYGLACTINE 910 BRAIDED

**BLUNT - 3/8 CIRCLE**

| NEEDLE                                                                                                 | EP  | USP | Length | Colour | Order Number |
|--------------------------------------------------------------------------------------------------------|-----|-----|--------|--------|--------------|
| <b>TERB 40 mm</b><br> | 3   | 2/0 | 90 cm  | violet | VC20TERB40   |
|                                                                                                        | 3.5 | 0   | 90 cm  | violet | VC00TERB40   |
|                                                                                                        | 4   | 1   | 90 cm  | violet | VC01TERB40   |
|                                                                                                        | 5   | 2   | 90 cm  | violet | VC02TERB40   |



ABSORBABLE SUTURES