




VENA NYLON



vena

Vena Medical Suppliers

www.sutures.com.eg

VENA NYLON

NONABSORBABLE SURGICAL SUTURE
NYLON MONOFILAMENT



Type

Monofilament

Composition

Extrusion of polyamide 6.0 or 6.6

Coating

None

Colour

Blue / Black

Tissue reaction

Minimal

Absorption

Non absorbable, gradually encapsulated by connective tissue.
The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).

BENEFITS

Excellent and permanent tissue support

Minimal tissue reaction

Soft passage through the tissues without sawing, tissue drag and trauma

*Absence of capillarity (no support of bacterial growth)
Excellent knot tie down*

*Good handling properties
Easy and secure knotting*

Guaranteed seal and product sterility

CHARACTERISTICS

High tensile strength, non absorbable

Synthetic

Monofilament with microscopic uniform diameter

Pliability and softness

Hermetically sealed packing

INDICATIONS

Daclon Blue

Skin closure / General surgery

Plastic surgery

Gastrointestinal surgery

Gynaecology


Obstetrics / Orthopaedics

Daclon Black

Plastic surgery / Ophthalmic surgery

VENA NYLON

NONABSORBABLE SURGICAL SUTURE
NYLON MONOFILAMENT



ROUND BODIED TAPER POINT - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHR 22 mm 	2	3/0	75 cm	blue	VN30OHR22
	3	2/0	75 cm	blue	VN20OHR22
OHR 26 mm 	1.5	4/0	75 cm	blue	VN40OHR26
	2	3/0	75 cm	blue	VN30OHR26
	3	2/0	75 cm	blue	VN20OHR26
	3.5	0	75 cm	blue	VN00OHR26
	4	1	75 cm	blue	VN01OHR26
OHR 30 mm 	3	2/0	75 cm	blue	VN20OHR30
	3.5	0	75 cm	blue	VN00OHR30
	4	1	75 cm	blue	VN01OHR30
OHR 40 mm 	3	2/0	75 cm	blue	VN20OHR40
	3.5	0	75 cm	blue	VN00OHR40
	4	1	75 cm	blue	VN01OHR40
	5	2	75 cm	blue	VN02OHR40


REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 19 mm 	1	5/0	75 cm	black	VN50TEX19
	1.5	4/0	45 cm	blue	VN40TEX19-45
	1.5	4/0	75 cm	blue	VN40TEX19
	2	3/0	45 cm	blue	VN30TEX19-45
	2	3/0	75 cm	blue	VN30TEX19
	3	2/0	75 cm	blue	VN20TEX19
TEX 24 mm 	1	5/0	75 cm	black	VN50TEX24
	1.5	4/0	45 cm	blue	VN40TEX24-45
	1.5	4/0	75 cm	blue	VN40TEX24
	2	3/0	45 cm	blue	VN30TEX24-45
	2	3/0	75 cm	blue	VN30TEX24
	3	2/0	75 cm	blue	VN20TEX24
	3.5	0	75 cm	blue	VN00TEX24
	4	1	75 cm	blue	VN01TEX24
TEX 26 mm 	2	3/0	45 cm	blue	VN30TEX26-45
	2	3/0	75 cm	blue	VN30TEX26
	3	2/0	45 cm	blue	VN20TEX26-45
	3	2/0	75 cm	blue	VN20TEX26





NONABSORBABLE SUTURES

VENA NYLON







NONABSORBABLE SURGICAL SUTURE
NYLON MONOFILAMENT



REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 30 mm 	2	3/0	75 cm	blue	VN30TEX30
	3	2/0	75 cm	blue	VN20TEX30
	3.5	0	75 cm	blue	VN00TEX30
	4	1	75 cm	blue	VN01TEX30
TEX 36 mm 	2	3/0	75 cm	blue	VN30TEX36
TEX 40 mm 	3	2/0	75 cm	blue	VN20TEX40
	3.5	0	75 cm	blue	VN00TEX40
	4	1	75 cm	blue	VN01TEX40
	5	2	90 cm	blue	VN02TEX40
TEX 45 mm 	3	2/0	100 cm	blue	VN20TEX45
	4	1	75 cm	blue	VN01TEX45
	5	2	90 cm	blue	VN02TEX45


REVERSE CUTTING - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHX 17 mm 	1.5	4/0	75 cm	blue	VN40OHX17
	2	3/0	75 cm	blue	VN30OHX17
OHX 22 mm 	2	3/0	75 cm	blue	VN30OHX22
OHX 26 mm 	2	3/0	75 cm	blue	VN30OHX26
	3	2/0	75 cm	blue	VN20OHX26
OHX 30 mm 	3.5	0	75 cm	blue	VN00OHX30
OHX 36 mm 	3.5	0	75 cm	blue	VN00OHX36
	4	1	75 cm	blue	VN01OHX36
	5	2	75 cm	blue	VN02OHX36
OHX 40 mm 	4	1	75 cm	blue	VN01OHX40




NONABSORBABLE SUTURES

VENA NYLON

NONABSORBABLE SURGICAL SUTURE
NYLON MONOFILAMENT



CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
 TEC 26 mm	2	3/0	75 cm	blue	VN30TEC26
 TEC 30 mm	3	2/0	75 cm	blue	VN20TEC30
 TEC 40 mm	3	2/0	75 cm	blue	VN20TEC40

NONABSORBABLE SUTURES