



VENA CRYL PDO



vena

Vena Medical Suppliers

www.sutures.com.eg

VENA CRYL PDO

ABSORBABLE SURGICAL SUTURE
POLYDIOXANONE MONOFILAMENT



Type	Monofilament
Composition	Polydioxanone, a polymer made from polyester poly (p-dioxanone)
Coating	None
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.

BENEFITS

- Outstanding strength retention for extended wound support*
- Minimal tissue reaction*
- Predictable absorption rate*
- Soft passage through the tissues without sawing, tissue drag and trauma*
- Absence of capillarity (no support of bacterial growth)*
- Excellent knot tie down*
- Good handling properties*
- Easy and secure knotting*
- Guaranteed seal and product sterility*

CHARACTERISTICS

- High and long-lasting tensile strength*
- Synthetic*
- Breaks down by hydrolysis*
- Monofilament with microscopic uniform diameter*
- Pliability and softness*
- Hermetically sealed packing*


INDICATIONS

- General soft tissue closing*
- General surgery*
- Gastrointestinal surgery*
- Ophthalmic surgery*
- Orthopaedics*
- Abdominal wall closure*
- Gynaecology*
- Plastic surgery*
- Urology*



VENA CRYL PDO

ABSORBABLE SURGICAL SUTURE
POLYDIOXANONE MONOFILAMENT



ROUND BODIED TAPER POINT - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHR 17 mm 	0.7	6/0	75 cm	violet	VD60OHR17
	1	5/0	75 cm	violet	VD50OHR17
	1.5	4/0	75 cm	violet	VD40OHR17
OHR 22 mm 	1.5	4/0	75 cm	violet	VD40OHR22
	2	3/0	75 cm	violet	VD30OHR22
OHR 26 mm 	1.5	4/0	75 cm	violet	VD40OHR26
	2	3/0	75 cm	violet	VD30OHR26
	3	2/0	75 cm	violet	VD20OHR26
	3.5	0	75 cm	violet	VD00OHR26
OHR 30 mm 	2	3/0	75 cm	violet	VD30OHR30
	3	2/0	75 cm	violet	VD20OHR30
	3.5	0	75 cm	violet	VD00OHR30
OHR 36 mm 	3	2/0	75 cm	violet	VD20OHR36
	3.5	0	90 cm	violet	VD00OHR36
	4	1	90 cm	violet	VD01OHR36
OHR 40 mm 	3.5	0	90 cm	violet	VD00OHR40
	4	1	90 cm	violet	VD01OHR40
	5	2	90 cm	violet	VD02OHR40

REVERSE CUTTING - 3/8 CIRCLE


NEEDLE	EP	USP	Length	Colour	Order Number
TEX 19 mm 	1	5/0	45 cm	violet	VD50TEX19
	1.5	4/0	75 cm	violet	VD40TEX19
	2	3/0	75 cm	violet	VD30TEX19
TEX 24 mm 	1.5	4/0	75 cm	violet	VD40TEX24
	2	3/0	75 cm	violet	VD30TEX24
	3	2/0	75 cm	violet	VD20TEX24
	3.5	0	75 cm	violet	VD00TEX24
TEX 26 mm 	2	3/0	75 cm	violet	VD30TEX26
	3	2/0	75 cm	violet	VD20TEX26

ABSORBABLE SUTURES



VENA CRYL PDO

ABSORBABLE SURGICAL SUTURE
POLYDIOXANONE MONOFILAMENT



REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 30 mm 	3	2/0	75 cm	violet	VD20TEX30
	3.5	0	75 cm	violet	VD00TEX30
TEX 40 mm 	4	1	90 cm	violet	VD01TEX40
	5	2	90 cm	violet	VD02TEX40

REVERSE CUTTING - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHX 36 mm 	3	2/0	75 cm	violet	VD20OHX36
	3.5	0	75 cm	violet	VD00OHX36
	4	1	75 cm	violet	VD01OHX36

ABSORBABLE SUTURES