



VENA PET



vena

Vena Medical Suppliers

www.sutures.com.eg

VENA PET

NONABSORBABLE SURGICAL SUTURE
POLYESTER BRAIDED (GREEN)



Type

Braided multifilament

Composition

Polyester - a polymer of polyethylene terephthalate

Coating

Silicone - Teflon

Colour

Green / White

Tissue reaction

Minimal

Absorption

Non absorbable, gradually encapsulated by connective tissue

BENEFITS

Excellent and permanent tissue support

Use in cardiovascular surgery (heart prosthesis)

Minimal tissue reaction

*Excellent handling properties
High flexibility
High tensile strength
Excellent knot security*

*Soft passage through the tissues with minimal sawing, tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action*

Guaranteed seal and product sterility

CHARACTERISTICS

Unsurpassed tensile strength, non absorbable

*Resistance to repeated bending
Synthetic*

Braided structure

Coated multifilament

Hermetically sealed packing

INDICATIONS

General soft tissue closing and/or ligation

Cardiovascular surgery





Ophthalmic surgery

General surgery / Skin closure

Orthopaedics / Gastrointestinal surgery

Obstetrics


ROUND BODIED TAPER POINT - 1/2 CIRCLE


NEEDLE	EP	USP	Length	Colour	Order Number
OHR 26 mm 	2	3/0	75 cm	green	VS30OHR26
	3	2/0	75 cm	green	VS20OHR26
	3.5	0	75 cm	green	VS00OHR26
OHR 30 mm 	3	2/0	75 cm	green	VS20OHR30
	3.5	0	75 cm	green	VS00OHR30
OHR 36 mm 	3	2/0	75 cm	green	VS20OHR36
	3.5	0	75 cm	green	VS00OHR36
	4	1	75 cm	green	VS01OHR36
OHR 40 mm 	4	1	75 cm	green	VS01OHR40
	5	2	75 cm	green	VS02OHR40



ROUND BODIED TAPER POINT - STRAIGHT

NEEDLE	EP	USP	Length	Colour	Order Number
STR 19 mm 	1.5	4/0	45 cm	green	VS40STR19
	2	3/0	50 cm	green	VS30STR19


REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 19 mm 	1	5/0	75 cm	green	VS50TEX19
	1.5	4/0	75 cm	green	VS40TEX19
	2	3/0	75 cm	green	VS30TEX19
	3	2/0	75 cm	green	VS20TEX19
TEX 24 mm 	2	3/0	75 cm	green	VS30TEX24
	3	2/0	75 cm	green	VS20TEX24
	3.5	0	75 cm	green	VS00TEX24


REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 26 mm 	3	2/0	75 cm	green	VS20TEX26
TEX 30 mm 	2	3/0	75 cm	green	VS30TEX30
	3	2/0	75 cm	green	VS20TEX30
	3.5	0	75 cm	green	VS00TEX30
	4	1	75 cm	green	VS01TEX30
TEX 36 mm 	3	2/0	75 cm	green	VS20TEX36
TEX 40 mm 	4	1	75 cm	green	VS01TEX40
	5	2	75 cm	green	VS02TEX40

REVERSE CUTTING - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHX 17 mm 	1.5	4/0	75 cm	green	VS40OHX17
	2	3/0	75 cm	green	VS30OHX17
OHX 22 mm 	2	3/0	75 cm	green	VS30OHX22