




VENA SILK



vena

Vena Medical Suppliers

www.sutures.com.eg

VENA SILK

NONABSORBABLE SURGICAL SUTURE
SILK BRAIDED



Type

Braided multifilament

Composition

Braided fibres from the cocoon of the silkworm

Coating

Waxes - Silicone

Colour

Black

Tissue reaction

Moderate

Absorption

Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.

BENEFITS

*Good and extended tissue support
for up to 3 month*

*Excellent handling properties
High flexibility
High tensile strength
Excellent knot security*

*Soft passage through the tissues
with minimal sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action*

Guaranteed seal and product sterility

CHARACTERISTICS

High tensile strength, non absorbable


Braided or twisted structure





Coated multifilament

Hermetically sealed packing


INDICATIONS

*General soft tissue closing
and/or ligation
General surgery
Skin closure / Ophthalmic surgery
Oral surgery / Gastrointestinal surgery
Gynaecology / Obstetrics*



ROUND BODIED TAPER POINT - 1/2 CIRCLE


NEEDLE	EP	USP	Length	Colour	Order Number
OHR 26 mm 	1	5/0	75 cm	black	VK50OHR26
	1.5	4/0	75 cm	black	VK40OHR26
	2	3/0	75 cm	black	VK30OHR26
	3	2/0	75 cm	black	VK20OHR26
	3.5	0	75 cm	black	VK00OHR26
	4	1	75 cm	black	VK01OHR26
OHR 30 mm 	2	3/0	75 cm	black	VK30OHR30
	3	2/0	75 cm	black	VK20OHR30
	3.5	0	75 cm	black	VK00OHR30
	4	1	75 cm	black	VK01OHR30
	5	2	75 cm	black	VK02OHR30
OHR 40 mm 	5	2	75 cm	black	VK02OHR40
OHR 45 mm 	4	1	75 cm	black	VK01OHR45
	5	2	75 cm	black	VK02OHR45




ROUND BODIED TAPER POINT - 3/8 CIRCLE


NEEDLE	EP	USP	Length	Colour	Order Number
TER 30 mm 	3	2/0	75 cm	black	VK20TER30





ROUND BODIED TAPER POINT - STRAIGHT

NEEDLE	EP	USP	Length	Colour	Order Number
STR 40 mm 	3	2/0	75 cm	black	VK20STR40


REVERSE CUTTING - 3/8 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
TEX 19 mm 	0.7	6/0	75 cm	black	VK60TEX19
	1	5/0	75 cm	black	VK50TEX19
	1.5	4/0	75 cm	black	VK40TEX19
	2	3/0	75 cm	black	VK30TEX19
	3	2/0	75 cm	black	VK20TEX19
	3.5	0	75 cm	black	VK00TEX19
	4	1	75 cm	black	VK01TEX19
TEX 24 mm 	1	5/0	75 cm	black	VK50TEX24
	1.5	4/0	75 cm	black	VK40TEX24
	2	3/0	75 cm	black	VK30TEX24
	3	2/0	75 cm	black	VK20TEX24
	3.5	0	75 cm	black	VK00TEX24
	4	1	75 cm	black	VK01TEX24
TEX 26 mm 	2	3/0	75 cm	black	VK30TEX26
	3	2/0	75 cm	black	VK20TEX26
TEX 30 mm 	1	5/0	75 cm	black	VK50TEX30
	1.5	4/0	75 cm	black	VK40TEX30
	2	3/0	75 cm	black	VK30TEX30
	3	2/0	75 cm	black	VK20TEX30
	3.5	0	75 cm	black	VK00TEX30
	4	1	75 cm	black	VK01TEX30
	5	2	75 cm	black	VK02TEX30
TEX 45 mm 	3	2/0	75 cm	black	VK20TEX45
	3.5	0	75 cm	black	VK00TEX45
	4	1	75 cm	black	VK01TEX45
	5	2	75 cm	black	VK02TEX45


REVERSE CUTTING - 1/2 CIRCLE

NEEDLE	EP	USP	Length	Colour	Order Number
OHX 17 mm 	1	5/0	75 cm	black	VK50OHX17
	1.5	4/0	75 cm	black	VK40OHX17
	2	3/0	75 cm	black	VK30OHX17
OHX 22 mm 	1	5/0	75 cm	black	VK50OHX22
	1.5	4/0	75 cm	black	VK40OHX22
	2	3/0	75 cm	black	VK30OHX22
	3	2/0	75 cm	black	VK20OHX22
OHX 26 mm 	1	5/0	75 cm	black	VK50OHX26
	1.5	4/0	75 cm	black	VK40OHX26
	2	3/0	75 cm	black	VK30OHX26
	3	2/0	75 cm	black	VK20OHX26
	3.5	0	75 cm	black	VK00OHX26
OHX 36 mm 	3.5	0	75 cm	black	VK00OHX36
	4	1	75 cm	black	VK01OHX36

CUTTING - STRAIGHT

NEEDLE	EP	USP	Length	Colour	Order Number
STC 40 mm 	3	2/0	75 cm	black	VK20STC40